

## CAUSTIC DIGESTER UNITS

ALKALINE HYDROLYSIS for  
DISPOSAL OF BIOLOGIC  
and INFECTIOUS WASTE

---

---

---

---

---

---

---

---

### ALKALINE HYDROLYSIS

- A natural process – animal carcasses buried in the earth are degraded by alkaline hydrolysis, expedited by the soil bacteria – a slow process
- Food in the intestine is digested to usable nutrients by alkaline hydrolysis, expedited by enzymes that operate at pH 7-8 at body temperature – a moderately fast process for relatively small amounts of tissue

---

---

---

---

---

---

---

---

### ALKALINE HYDROLYSIS

Uses strong alkali (pH 14) to solubilize and hydrolyze tissue, expedited by heat in a pressurized vessel – **a very fast process for large amounts of tissue**

Digests and sterilizes in one operation

Generates an EPA neutral solution of amino acids, peptides, sugars, soaps, and electrolytes that is suitable for release to a sanitary sewer or for use as fertilizer or as feedstock for biogas or biodiesel production

---

---

---

---

---

---

---

---

## ALKALINE HYDROLYSIS

**Proteins:** Hydrolyzed to amino acids and small peptides by breaking of amide (peptide) bonds; carbohydrates clipped from glycoproteins; specific amino acids destroyed

**Fats:** Ester bonds between fatty acids and glycerol hydrolyzed, yielding soaps; glycolipids and polyunsaturated fats also "destroyed"

**Nucleic acids:** Phosphodiester bonds of nucleotide chains hydrolyzed, RNA rapidly, DNA more slowly

---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

---

---

## CAUSTIC DIGESTION

- Can reduce volume and weight of biologic waste by as much as 97%
- Treats fresh, frozen, and fixed tissue equally well
- Sterilizes any non-biologic materials present in waste

---

---

---

---

---

---

---

---

---

---

## The result of digestion:

330 g rat and its littermate (in dish) after caustic digestion




---

---

---

---

---

---

---

---

---

---

## CAUSTIC DIGESTION

- Destroys all pathogens, *including prions*
- Converts fixatives, cytotoxic agents, and other toxins to harmless, biodegradable derivatives
- Releases radionuclides from tissue into an aqueous solution suitable for release to a sanitary sewer under 10CFR20

---

---

---

---

---

---

---

---

## Efficacy Testing Results

Biological Indicator	Before Digestion	After Digestion
<i>Aspergillus fumigatus</i>	9.0 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Bacillus subtilis</i> (vegetative)	6.0 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Pseudomonas aeruginosa</i>	2.3 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Giardia</i> cysts	2.0 x 10 <sup>6</sup>	zero per ml examined
<i>Mycobacterium bovis</i> BCG	4.0 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Giardia muris</i>	89% excystation	cysts completely destroyed
MS-2 bacteriophage	1.0 x 10 <sup>6</sup> PFU/ml	zero per ml tested
<i>Staphylococcus aureus</i>	4.0 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Mycobacterium fortuitum</i>	6.0 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Candida albicans</i>	1.4 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Mycobacterium terrae</i>	3 x 10 <sup>6</sup> CFU/ml	zero per ml tested
<i>Bacillus stearothermophilus</i>	3M spore strips	no spores detected

---

---

---

---

---

---

---

---

**All Select Agents as well as protein and non-protein toxins, are destroyed by alkaline hydrolysis at elevated temperature (i.e., Caustic Digestion)**

---

---

---

---

---

---

---

---

## CAUSTIC DIGESTERS

**Prototype** -1993 Albany Med

**Development** - 1994-2006 WR<sup>2</sup>

**Re-engineering and improvement** - 2006--- PRI-BIO

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---




---

---

---

---

---

---

---

---

Early Design for a 3500 lb Capacity Unit



---

---

---

---

---

---

---

---

Loading partial bovine carcass into the Tissue Digester



---

---

---

---

---

---

---

---

30 lb  
Capacity  
for  
Small  
Laboratory



---

---

---

---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

## PRI-BIO CDU

- Engineering Improvements:
  - Industrial Grade Components and Parts
  - Direct Steam Injection Heating and other options
  - High-Temp Magnetically-Driven Pumps
  - Hydraulic Ring Closures Mounted on the Vessel
  - Full Insulation of External Piping
  - External Heat Exchanger Cooling or other options
  - Automated pH and Odor Reduction
  - Automated Batch Reporting and Record Keeping
  - Effluent BOD Reduction Option
  - *All Manufacturing Done at PRI Plant*

---

---

---

---

---

---

---

---





Paired 500 lb Capacity CDU for NIH

---

---

---

---

---

---

---

---



500 lb Capacity CDU for NIH

---

---

---

---

---

---

---

---

5,000 lb  
Capacity  
CDU for  
Cornell  
CVM  
in PRI  
Plant



---

---

---

---

---

---

---

---



---

---

---

---

---

---

---

---



---

---

---

---

---

---

---

---



---

---

---

---

---

---

---

---



Sliding Skid Into Place on  
PRI-Engineered Rail System

---

---

---

---

---

---

---

---



Lowering Vessel Onto Skid



---

---

---

---


---

---

---

---

Digester/Sterilizer Combo Pass-Through



---

---

---

---

---

---

---

---

80 lb Capacity  
CDU for Tufts  
University  
Manual Latching  
Heavy Flat Lid



---

---

---

---

---

---

---

---

Tufts CDU  
With Basket  
Lifted



---

---

---

---

---

---

---

---

Tufts CDU  
With  
Safety  
Insulation  
Jacket Over  
Lid



---

---

---

---

---

---

---

---



All External Piping is Safety Insulated

---

---

---

---

---

---

---

---

PRI-BIO CDU are matched to the needs of the customer in sizes ranging from 80 lb to 10,000 lb capacity.

---

---

---

---

---

---

---

---